

Product Datasheet

USP28 Blocking Peptide (orb706064)

Description	The peptide is used to block Anti-USP28 antibody < a href="https://www.biorbyt.com/usp28-antibody-or..."
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	USP28
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q96RU2 , Q5I043
Entrez	235323 , 57646
Expiration Date	12 months from date of receipt.