

Product Datasheet

RGPD1/2/3/4/5/8 Blocking Peptide (orb706125)

Description	The peptide is used to block Anti-RGPD1/2/3/4/5/8 antibody < a href="https://www.biorbyt.com/rgpd1-a...
Reactivity	Human, Rat
Tested Applications	Blocking
Target	RGPD1; RGPD2; RGPD3; RGPD4; RGPD5; RGPD8
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P0DJD0 , P0DJD1 , O14715 , Q99666 , A6NKT7 , Q7Z3J3
Entrez	84220 , 729540 , 729857
Expiration Date	12 months from date of receipt.