
Product Datasheet

CaMK2 alpha/beta/delta (phospho-T305) Blocking Peptide (orb706245)

Description

The peptide is used to block Anti-

antibody < a
href="https://www.bior...

Reactivity

Human, Mouse, Porcine, Rabbit, Rat, Zebrafish

Tested Applications

Blocking

Target

CAMK2A

Form/Appearance

Lyophilized powder

Storage

Shipped at 4°C. Store at -20°C for one year.

Note

For research use only

Application notes

Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.

Purity

> 85%

Source

Synthetic

Uniprot ID

[Q13557](#), [P15791](#), [P08413](#), [P11798](#), [Q9UQM7](#), [P28652](#), [P11275](#), [Q6PHZ2](#), [Q13554](#)

Entrez

[108058](#), [12322](#), [25400](#), [816](#), [24246](#), [12323](#), [24245](#), [817](#)

Expiration Date

12 months from date of receipt.