

Product Datasheet

C/EBP alpha (phospho-S21) Blocking Peptide (orb706248)

Description	The peptide is used to block Anti-C/EBP alpha (pS21) antibody < a href="https://www.biorbyt.com/cebp..."
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CEBPA
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P05554 , P53566 , P49715
Entrez	1050 , 24252 , 12606
Expiration Date	12 months from date of receipt.