

Product Datasheet

C/EBP epsilon (phospho-T74) Blocking Peptide (orb706258)

| | |
|----------------------------|---|
| Description | The peptide is used to block Anti-C/EBP epsilon (pT74) antibody < a href="https://www.biorbyt.com/ce... |
| Reactivity | Human, Mouse, Rat, Sheep |
| Tested Applications | Blocking |
| Target | CEBPE |
| Form/Appearance | Lyophilized powder |
| Storage | Shipped at 4°C. Store at -20°C for one year. |
| Note | For research use only |
| Application notes | Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours. |
| Purity | > 85% |
| Source | Synthetic |
| Uniprot ID | P56261 , Q15744 , Q6PZD9 |
| Entrez | 25410 , 110794 , 1053 |
| Expiration Date | 12 months from date of receipt. |