

Product Datasheet

Ubiquitin (phospho-S65) Blocking Peptide (orb706272)

Description	The peptide is used to block Anti-Ubiquitin (pS65) antibody < a href="https://www.biorbyt.com/uba52-...>
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	UBA52; RPS27A; UBB; UBC
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P62987 , P62979 , P0CG47 , P0CG48
Entrez	7314 , 6233 , 7311 , 7316
Expiration Date	12 months from date of receipt.