

Product Datasheet

MARCKS (phospho-S158) Blocking Peptide (orb706276)

Description	The peptide is used to block Anti-MARCKS (pS158) antibody < a href="https://www.biorbyt.com/marcks-p..."
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	MARCKS
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P26645 , P29966 , P30009
Entrez	17118 , 4082
Expiration Date	12 months from date of receipt.