

Product Datasheet

DSCR1 (phospho-S108) Blocking Peptide (orb706287)

Description	The peptide is used to block Anti-DSCR1 (pS108) antibody < a href="https://www.biorbyt.com/rcan1-pho...">
Reactivity	Human, Mouse, Rat, Sheep
Tested Applications	Blocking
Target	RCAN1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9JHG6 , P53805 , Q6IN33
Entrez	1827 , 54720 , 266766
Expiration Date	12 months from date of receipt.