

Product Datasheet

AMPK gamma 1/2/3 Blocking Peptide (orb706324)

Description	The peptide is used to block Anti-AMPK gamma 1/2/3 antibody < a href="https://www.biorbyt.com/prkag1..."
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	PRKAG1; PRKAG2; PRKAG3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q91WG5 , O54950 , Q9UGI9 , P80385 , Q9UGJ0 , P54619 , Q8BGM7
Entrez	241113 , 53632 , 25520 , 108099 , 5571 , 19082 , 51422
Expiration Date	12 months from date of receipt.