

Product Datasheet

MKK5 (phospho-S311) Blocking Peptide (orb706334)

Description	The peptide is used to block Anti-MKK5 (pS311) antibody < a href="https://www.biorbyt.com/map2k5-pho..."
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	MAP2K5
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9WVS7 , Q62862 , Q13163
Entrez	29568 , 23938 , 5607
Expiration Date	12 months from date of receipt.