
Product Datasheet

Histone H2B (AcK85) Blocking Peptide (orb706359)

Description

The peptide is used to block Anti-Histone

[href=https://www.biorbyt.com/hist...](https://www.biorbyt.com/hist...)

Reactivity

Human, Mouse, Rat

Tested Applications

Blocking

Target

HIST1H2BC; HIST1H2BE; HIST1H2BF;
HIST1H2BG; HIST1H2BI; HIST1H2BD;
HIST1H2BH; HIST1H2BK; HIST1H2BL;
HIST1H2BM; HIST1H2BN; HIST2H2BF;
H2BFS

Form/Appearance

Lyophilized powder

Storage

Shipped at 4°C. Store at -20°C for one year.

Note

For research use only

Application notes

Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.

Purity

> 85%

Source

Synthetic

Uniprot ID

[Q99880](#), [Q99877](#), [Q99879](#), [Q93079](#),
[O60814](#), [P62807](#), [Q6ZWY9](#), [P58876](#),
[Q64475](#), [Q64478](#), [P23527](#), [P33778](#),
[P06899](#)

Entrez

[8344](#), [8341](#), [319182](#), [440689](#), [319179](#),
[319178](#), [85236](#), [8342](#), [8340](#), [8339](#), [8345](#)
, [8343](#), [8347](#), [8346](#), [3017](#)

Expiration Date

12 months from date of receipt.