

Product Datasheet

MYOD1 (phospho-S200) Blocking Peptide (orb706368)

Description	The peptide is used to block Anti-MYOD1 (pS200) antibody < a href=https://www.biorbyt.com/myod1-ps20...>
Reactivity	Human, Mouse, Porcine, Rat, Sheep, Zebrafish
Tested Applications	Blocking
Target	MYOD1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P15172 , Q02346 , P10085
Entrez	337868 , 4654 , 17927
Expiration Date	12 months from date of receipt.