

Product Datasheet

NF-kappaB p105 (phospho-S903) Blocking Peptide (orb706370)

| | |
|----------------------------|---|
| Description | The peptide is used to block Anti-NF-kappaB p105 (pS903) antibody < a href=https://www.biorbyt.com/n...> |
| Reactivity | Human, Rat |
| Tested Applications | Blocking |
| Target | NFKB1 |
| Form/Appearance | Lyophilized powder |
| Storage | Shipped at 4°C. Store at -20°C for one year. |
| Note | For research use only |
| Application notes | Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours. |
| Purity | > 85% |
| Source | Synthetic |
| Uniprot ID | P19838, Q63369 |
| Entrez | 4790 |
| Expiration Date | 12 months from date of receipt. |