

Product Datasheet

VEGFR2 (phospho-Y1059) Blocking Peptide (orb706374)

Description	The peptide is used to block Anti-VEGFR2 (pY1059) antibody < a href=https://www.biorbyt.com/vegfr2-p...
Reactivity	Human, Mouse, Zebrafish
Tested Applications	Blocking
Target	KDR
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P35918, P35968
Entrez	3791
Expiration Date	12 months from date of receipt.