

Product Datasheet

Cyclin D1 (phospho-S90) Blocking Peptide (orb706386)

Description	The peptide is used to block Anti-Cyclin D1 (pS90) antibody < a href=https://www.biorbyt.com/cyclin-...
Reactivity	Human, Mouse, Rat, Sheep
Tested Applications	Blocking
Target	CCND1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P25322 , P39948 , P24385
Entrez	12443 , 595 , 58919
Expiration Date	12 months from date of receipt.