

Product Datasheet

Focal Adhesion Kinase (phospho-S843) Blocking Peptide (orb706407)

Description	The peptide is used to block Anti-Focal Adhesion Kinase (pS843) antibody < a href=https://www.biorby...
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	PTK2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P34152 , O35346 , Q05397
Entrez	5747 , 14083 , 25614
Expiration Date	12 months from date of receipt.