

Product Datasheet

MERTK/TYRO3 Blocking Peptide (orb706429)

Description	The peptide is used to block Anti-MERTK/TYRO3 antibody < a href=https://www.biorbyt.com/mertk-tyro3-...
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	MERTK; TYRO3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P55146 , Q60805 , Q06418 , P55144 , Q12866 , P57097
Entrez	10461 , 7301 , 17289 , 25232 , 22174 , 65037
Expiration Date	12 months from date of receipt.