

Product Datasheet

CD66a/b/c Blocking Peptide (orb706532)

Description	The peptide is used to block Anti-CD66a/b/c antibody < a href=https://www.biorbyt.com/cd66a-b-c-anti...>
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	CEACAM1; CEACAM6; CEACAM8
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P40199 , P13688 , P31997
Entrez	1088 , 4680 , 634
Expiration Date	12 months from date of receipt.