

Product Datasheet

Synaptotagmin (phospho-S309/306) Blocking Peptide (orb706538)

Description	The peptide is used to block Anti-Synaptotagmin (pS309/306) antibody < a href=https://www.biorbyt.co...
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	SYT1; SYT2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P21579 , P29101 , P21707 , P46096 , P46097 , Q8N9I0
Entrez	20979 , 20980 , 127833 , 24805 , 25716 , 6857
Expiration Date	12 months from date of receipt.