

Product Datasheet

ABL1/2 (phospho-Y393/439) Blocking Peptide (orb706539)

Description	The peptide is used to block Anti-ABL1/2 (pY393/439) antibody < a href=https://www.biorbyt.com/abl1-...
Reactivity	Human, Mouse, Primate
Tested Applications	Blocking
Target	ABL1; ABL2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P00519, P42684, P00520, Q4JIM5
Entrez	11350, 27, 25
Expiration Date	12 months from date of receipt.