

Product Datasheet

FOXO1/3/4 (phospho-T24/32) Blocking Peptide (orb706541)

Description	The peptide is used to block Anti-FOXO1/3/4 (pT24/32) antibody < a href=https://www.biorbyt.com/foxo...
Reactivity	Human, Mouse, Porcine, Rat, Zebrafish
Tested Applications	Blocking
Target	FOXO1; FOXO3; FOXO4
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9R1E0 , Q9WVH3 , G3V7R4 , Q9WVH4 , O43524 , Q12778 , P98177
Entrez	54601 , 56458 , 2308 , 56484 , 4303 , 84482
Expiration Date	12 months from date of receipt.