

Product Datasheet

c-RAF (phospho-Y341) Blocking Peptide (orb706542)

Description	The peptide is used to block Anti-c-RAF (pY341) antibody < a href=https://www.biorbyt.com/c-raf-py34...>
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	RAF1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q99N57 , P04049 , P11345
Entrez	24703 , 5894 , 110157
Expiration Date	12 months from date of receipt.