

Product Datasheet

p23 (phospho-S113) Blocking Peptide (orb706543)

Description	The peptide is used to block Anti-p23 (pS113) antibody < a href=https://www.biorbyt.com/p23-ps113-an... an...
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	PTGES3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9R0Q7 , Q15185 , P83868
Entrez	102641464 , 362809 , 10728 , 56351
Expiration Date	12 months from date of receipt.