

Product Datasheet

RAC1/2/3 (phospho-S71) Blocking Peptide (orb706550)

Description	The peptide is used to block Anti-RAC1/2/3 (pS71) antibody < a href=https://www.biorbyt.com/rac1-2-3...>
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	RAC1; RAC2; RAC3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8CFN2 , P60764 , P60766 , P63001 , P60763 , Q6RUV5 , P63000 , Q05144 , P15153
Entrez	64465 , 998 , 12540 , 5881 , 19354 , 170758 , 363875 , 5879 , 19353 , 5880
Expiration Date	12 months from date of receipt.