

Product Datasheet

HSP27 (phospho-S78) Blocking Peptide (orb706551)

Description	The peptide is used to block Anti-HSP27 (pS78) antibody < a href=https://www.biorbyt.com/hsp27-ps78-...>
Reactivity	Human, Rat
Tested Applications	Blocking
Target	HSPB1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P04792 , P42930
Entrez	24471 , 3315
Expiration Date	12 months from date of receipt.