

Product Datasheet

HSF1 (phospho-S121) Blocking Peptide (orb706555)

| | |
|----------------------------|---|
| Description | The peptide is used to block Anti-HSF1 (pS121) Antibody (orb665666) reactivity. |
| Reactivity | Human, Mouse, Rat |
| Tested Applications | Blocking |
| Target | HSF1 |
| Form/Appearance | Lyophilized powder |
| Storage | Shipped at 4°C. Store at -20°C for one year. |
| Note | For research use only |
| Application notes | Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours. |
| Purity | > 85% |
| Source | Synthetic |
| Uniprot ID | P38532, Q00613 |
| Entrez | 15499, 3297 |
| Expiration Date | 12 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)