

Product Datasheet

BLNK (phospho-Y84) Blocking Peptide (orb706559)

Description	The peptide is used to block Anti-BLNK (pY84) antibody < a href=https://www.biorbyt.com/blnk-py84-an...>
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	BLNK
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q4KM52 , Q9QUN3 , Q8WV28
Entrez	499356 , 29760 , 17060
Expiration Date	12 months from date of receipt.