

Product Datasheet

GAB2 (phospho-Y643) Blocking Peptide (orb706564)

Description	The peptide is used to block Anti-GAB2 (pY643) antibody < a href=https://www.biorbyt.com/gab2-py643-...>
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	GAB2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9EQH1 , Q9UQC2 , Q9Z1S8
Entrez	84477 , 9846
Expiration Date	12 months from date of receipt.