

## Product Datasheet

### GABBR2 (phospho-S893) Blocking Peptide (orb706576)

|                            |                                                                                                                                                                                                                                                                                                                 |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | The peptide is used to block Anti-GABBR2 (pS893) antibody < a href=https://www.biorbyt.com/gabbr2-ps...                                                                                                                                                                                                         |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                                                                                                                                                                                                               |
| <b>Tested Applications</b> | Blocking                                                                                                                                                                                                                                                                                                        |
| <b>Target</b>              | GABBR2                                                                                                                                                                                                                                                                                                          |
| <b>Form/Appearance</b>     | Lyophilized powder                                                                                                                                                                                                                                                                                              |
| <b>Storage</b>             | Shipped at 4°C. Store at -20°C for one year.                                                                                                                                                                                                                                                                    |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                                                                                           |
| <b>Application notes</b>   | Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours. |
| <b>Purity</b>              | > 85%                                                                                                                                                                                                                                                                                                           |
| <b>Source</b>              | Synthetic                                                                                                                                                                                                                                                                                                       |
| <b>Uniprot ID</b>          | <b>Q80T41, O88871, O75899</b>                                                                                                                                                                                                                                                                                   |
| <b>Entrez</b>              | <b>83633, 9568, 242425</b>                                                                                                                                                                                                                                                                                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                                                                                                                                 |