

Product Datasheet

TRK B (phospho-Y817) Blocking Peptide (orb706578)

Description	The peptide is used to block Anti-TRK B (pY817) antibody < a href=https://www.biorbyt.com/trk-b-py81...>
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	NTRK2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q63604 , Q16620 , P15209
Entrez	18212 , 25054 , 4915
Expiration Date	12 months from date of receipt.