

Product Datasheet

FGFR4 (phospho-Y642) Blocking Peptide (orb706581)

Description	The peptide is used to block Anti-FGFR4 (pY642) antibody < a href=https://www.biorbyt.com/fgfr4-py64...
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	FGFR4
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q498D6 , P22455 , Q03142
Entrez	25114 , 14186 , 2264
Expiration Date	12 months from date of receipt.