

## Product Datasheet

# Cytochrome P450 2C8/9/18/19 Blocking Peptide (orb706651)

<b>Description</b>	The peptide is used to block Anti-Cytochrome P450 2C8/9/18/19 antibody < a href=https://www.biorbyt....
<b>Reactivity</b>	Human, Mouse, Primate, Rat
<b>Tested Applications</b>	Blocking
<b>Target</b>	CYP2C8; CYP2C9; CYP2C18; CYP2C19
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year.
<b>Note</b>	For research use only
<b>Application notes</b>	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
<b>Purity</b>	> 85%
<b>Source</b>	Synthetic
<b>Uniprot ID</b>	<b>P10632, P33261, P11712, P33260</b>
<b>Entrez</b>	<b>1558, 1562, 1557, 1559</b>
<b>Expiration Date</b>	12 months from date of receipt.