

Product Datasheet

EPHB1/2/3 Blocking Peptide (orb706750)

Description	The peptide is used to block Anti-EPHB1/2/3 antibody < a href=https://www.biorbyt.com/ephb1-2-3-anti...>
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	EPHB1; EPHB2; EPHB3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P54763 , P54762 , Q8CBF3 , P54753 , P29323 , P54754
Entrez	13845 , 270190 , 2049 , 2047 , 2048
Expiration Date	12 months from date of receipt.