

---

## Product Datasheet

### Human BCL2 protein (orb706811)

**Description**

Recombinant Human BCL2 Protein is

target gene encoding ...

**Reactivity**

Human

**Endotoxins**

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

**Tested Applications**

ELISA, MS, SDS-PAGE, WB

**Target**

BCL2

**Preservatives**

Liquid in a 0.2 μm filtered solution of 20mM HEPES, 150mM NaCl, 10%Glycerol, pH 8.0.

**Form/Appearance**

Liquid in a 0.2 μm filtered solution of 20mM HEPES, 150mM NaCl, 10%Glycerol, pH 8.0.

**Storage**

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

**Note**

For research use only

**Application notes**

Please minimize freeze-thaw cycles.

**Protein Sequence**

 MAHAGRTGYD NREIVMKYIH YKLSQRGYEW  
 DAGDVGAAPP GAAPAGIFS SQPGHTPPHA  
 ASRDPVARTS PLQTPAAPGA AAGPALSPVP  
 PVVHLTLRQA GDDFSRRYRR DFAEMSSQLH  
 LTPFTARGRF ATVVEELFRD GVNWGRIVAF  
 FEEGGVMCVE SVNREMSPLV DNIALWMTEY  
 LNRHLHTWIQ DNGGWDAFVE LYGPMSMRPLF  
 DHHHHHH

**Purity**

Greater than 95% as determined by reducing SDS-PAGE.

**Source**

E. coli

**Uniprot ID**
**P10415**
**Entrez**
**596**
**Expiration Date**

6 months from date of receipt.