

## Product Datasheet

### Human CD38 protein (orb706884)

**Description**

Recombinant Human CD38 Protein is

and the target gene encodi...

**Reactivity**

Human

**Endotoxins**

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

**Tested Applications**

ELISA, MS, SDS-PAGE, WB

**Target**

CD38

**Preservatives**

Liquid in a 0.2 μm filtered solution of PBS, pH 7.4.

**Form/Appearance**

Liquid in a 0.2 μm filtered solution of PBS, pH 7.4.

**Storage**

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

**Note**

For research use only

**Application notes**

Please minimize freeze-thaw cycles.

**Protein Sequence**

VPRWRQQWSG PGTTKRFPET VLARCVKYTE  
 IHPEMRHVDC QSVWDAFKGA FISKHPCNIT  
 EEDYQPLMKL GTQTVPCNKI LLWSRIKDLA  
 HQFTQVQRDM FTLEDTLG Y LADDLTWCGE  
 FNTSKINYQS CPDWRKDCSN NPVSFVKTV  
 SRRFAEAACD VVHVMLNGSR SKIFDKNSTF  
 GSVEVHNLQP EKVQTLEAWV IHGGREDSRD  
 LCQDPTIKEL ESIISKRNIQ FSCKNIYRPD  
 KFLQCVKNPE DSSCTSEIHH HHHH

**Purity**

Greater than 95% as determined by reducing SDS-PAGE.

**Source**

Human cells

**Uniprot ID**
**P28907**
**Entrez**
**952**
**Expiration Date**

6 months from date of receipt.