
Product Datasheet

Human MET protein (orb706910)

Description

Recombinant Human c-Met Protein is

and the target gene encod...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

MET

Preservatives

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

 Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

ECKEALAKSE MNVNMKYQLP NFTAETPIQN
 VILHEHHIFL GATNYIYVLN EEDLQKVAEY
 KTGPVLEHPD CFPCQDCSSK ANLSGGVWKD
 NINMALVVDY YDDQLISCG SVNRGTCQRH
 VFPHNHTADI QSEVHCIFSP QIEEPSQCPD
 CVVSALGAKV LSSVKDRFIN FVVGNTINSS
 YFPDHPLHSI SVRRLKETKD GFMFLTDQSY
 IDVLPFRDS YPIKYVHAPE SNNFIYFLTV
 QRETLDAQTF HTRIIRFCSI NSGLHXYMEM
 PLECILTEKR KKRSTKKEVF NILQAAYVSK
 PGAQLARQIG ASLNDLILFG VFAQSKPDSA
 EPMDRSAMCA FPIKYVNDFF NKIVNKNVVR
 CLQHFYGNH EHCNRTLLR NSSGCEARRD
 EYRTEFTAL QRVDLFMGQF SEVLLTSIST
 FIKGDLTIAN LGTSEGRFMQ VVVSRSRGPST
 PHVNFLLDLH PVSPEVIVEH TLNQNGYTLV
 ITGKKITKIP LNGLGHHHHH H