
Product Datasheet

Human FGF2 protein (orb706938)

Description

Recombinant Human FGF2 Protein is

target gene encoding ...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

FGF2

Preservatives

Lyophilized from a 0.2 μm filtered solution of 20mM Tris, 150mM NaCl, pH 7.5.

Form/Appearance

Lyophilized from a 0.2 μm filtered solution of 20mM Tris, 150mM NaCl, pH 7.5.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

 Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

 MAAGSITLTP ALPEDGGSGA FPPGHFKDPK
 RLYCKNGGFF LRIHPDGRVD GVREKSDPHI
 KLQLQAEERG VVSIKGVCAN RYLAMKEDGR
 LLASKCVTDE CFFFERLESN NYNTYRSRKY
 TSWYVALNRT GQYKLGSKTG PGQKAILFLP
 MSAKS

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

E. coli

Uniprot ID
[P09038](#)
Entrez
[2247](#)
Expiration Date

6 months from date of receipt.