

---

## Product Datasheet

### Human HBEGF protein (orb706952)

**Description**

Recombinant Human HBEGF Protein is

and the target gene encod...

**Reactivity**

Human

**Endotoxins**

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

**Tested Applications**

ELISA, MS, SDS-PAGE, WB

**Target**

HBEGF

**Preservatives**

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

**Form/Appearance**

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

**Storage**

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

**Note**

For research use only

**Application notes**

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in ddH4O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

**Protein Sequence**

LVTGESLERL RRGLAAGTSN PDPPTVSTDQ  
 LLPLGGGRDR KVRDLQEADL DLLRVTLSSK  
 PQALATPNKE EHGKRKKKGK GLGKKRDPCL  
 RKYKDFCIHG ECKYVKELRA PSCICHPGYH  
 GERCHGLSLH HHHHH

**Purity**

Greater than 95% as determined by reducing SDS-PAGE.

**Source**

Human cells

**Uniprot ID**
[Q99075](#)
**Entrez**
[1839](#)
**Expiration Date**

6 months from date of receipt.