

---

## Product Datasheet

### Human IL1 Receptor protein (orb706974)

**Description**

Recombinant Human IL-1R1 Protein is

system and the target gene enco...

**Reactivity**

Human

**Endotoxins**

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

**Tested Applications**

ELISA, MS, SDS-PAGE, WB

**Target**

IL1R1

**Preservatives**

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

**Form/Appearance**

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

**Storage**

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

**Note**

For research use only

**Application notes**

 Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in ddH<sub>2</sub>O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

**Protein Sequence**

```

LEADKCKERE EKILVSSAN EIDVRPCPLN
PNEHKGITIW YKDDSKTPVS TEQASRIHQH
KEKLWFVPAK VEDSGHYCV VRNSSYCLRI
KISAKFVENE PNLCYNAQAI FKQKLPVAGD
GGLVCPYMEF FKNENNELPK LQWYKDCCKPL
LLDNIHFSGV KDRLIVMNVA EKHRGNYTCH
ASYTYLGKQY PITRVIEFIT LEENKPTRPV
IVSPANETME VDLGSQIQLI CNVTGQLSDI
AYWKWNGSVI DEDDPVLGED YYSVENPANK
RRSTLITVLN ISEIESRFYK HPFTCFAKNT
HGIDAAYIQL IYPVTIEGRD MDPKSCDKTH
TCPPCPAPEL LGGPSVFLFP PKPKDTLMIS
RTPEVTCVVV DVSHEDPEVK FNWYVDGVEV
HNAKTKPREE QYNSTYRVVS VLTVLHQDWL
NGKEYKCKVS NKALPAPIEK TISKAKGQPR
EPQVYTLPPS RDELTKNQVS LTCLVKGFYP
SDIAVEWESN GQPENNYKTT PPVLDSDGSF
ELYSKI TVDK SRWOOGNVFS CSVMHFAI HN
  
```

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

 Email: [info@biorbyt.com](mailto:info@biorbyt.com) | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive&lt;br&gt;Research Triangle Park&lt;br&gt;Durham, North Carolina&lt;br&gt;27709. United States

 Email: [info@biorbyt.com](mailto:info@biorbyt.com) | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558