
Product Datasheet

Human Mesothelin protein (orb707004)

Description

Recombinant Human Mesothelin Protein is

system and the target gene ...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

MSLN

Preservatives

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

LAGETGQEEA PLDGVLANPP NISSLSRQL
 LGFPCAEVSG LSTERVRELA VALAQKNVKL
 STEQLRCLAH RLSEPPEDLD ALPLDLLFL
 NPDAFSGPQA CTRFFSRITK ANVDLLPRGA
 PERQRLLPAA LACWGVGRSL LSEADVRLG
 GLACDLPGRF VAESAEVLLP RLVSCPGPLD
 QDQEEAARAA LQGGGPPYGP PSTWSVSTMD
 ALRGLLPVLG QPIIRSIPQG IVAAWRQRSS
 RDPSWRQPER VDDIEGRMDE PKSCDKTHTC
 PPCPAPELLG GPSVFLFPPK PKDTLMISRT
 PEVTCVVVDV SHEDPEVKFN WYVDGVEVHN
 AKTKPREEQY NSTYRVVSVL TVLHQDWLNG
 KEYKCKVSNK ALPAPIEKTI SKAKGQPREP
 QVYTLPPSRE EMTKNQVSLT CLVKGFYPSD
 IAVEWESNGQ PENNYKTPPP VLDSGDGFFL
 YSKLTVDKSR WQQGNVFSCS VMHEALHNNHY
 TQKSLSLSPG K