

---

## Product Datasheet

### Human VSIG8 protein (orb707074)

**Description**

Recombinant Human VSIG8 Protein is

and the target gene encod...

**Reactivity**

Human

**Endotoxins**

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

**Tested Applications**

ELISA, MS, SDS-PAGE, WB

**Target**

VSIG8

**Preservatives**

Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

**Form/Appearance**

Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

**Storage**

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

**Note**

For research use only

**Application notes**

 Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in ddH<sub>2</sub>O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

**Protein Sequence**

 VRINGDGQEV LYLAEGDNVR LGCPYVLDPE  
 DYGPNGLDIE WMQVNSDPAH HRENVFLSYQ  
 DKRINHGSLP HLQQRVRF AA SDPSQYDASI  
 NLMNLQVSDT ATYECRVKKT TMATRKVIVT  
 VQARPAVPMC WTEGHMTYGN DVVLKCYASG  
 GSQPLSYKWA KISGHHYPYR AGSYTSQHSY  
 HSELSYQESF HSSINQGLNN GDLVLKDISR  
 ADDGLYQCTV ANNVGYSVCV VEVKVSRSRR  
 IGVDDHHHHHH

**Purity**

Greater than 95% as determined by reducing SDS-PAGE.

**Source**

Human cells

**Uniprot ID**
[Q5VU13](#)
**Entrez**
[391123, 284677](#)