

## Product Datasheet

### Mouse CSF1 protein (orb707253)

**Description**

Recombinant Mouse M-CSF Protein is

and the target gene encod...

**Reactivity**

Mouse

**Endotoxins**

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

**Tested Applications**

ELISA, MS, SDS-PAGE, WB

**Target**

CSF1

**Preservatives**

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

**Form/Appearance**

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

**Storage**

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

**Note**

For research use only

**Application notes**

 Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in ddH<sub>2</sub>O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

**Protein Sequence**

 KEVSEHCSHM IGNGHLKVLQ QLIDSQMETS  
 CQIAFEFVDQ EQLDDPVCYL KKAFFLVQDI  
 IDETMRFKDN TPNANATERL QELSNNLNSC  
 FTKDYEEQNK ACVRTFHETP LQLEKIKNF  
 FNETKNLLEK DWNIFTKNCN NSFACSSRD  
 VVTKPDCNCL YPKATPSSDP ASASPHQPPA  
 PSMAPLAGLA WDDSQRTEGS SLLPSELPLR  
 IEDPGSAKQR PPRSTCQMLE VDHHHHHH

**Purity**

Greater than 95% as determined by reducing SDS-PAGE.

**Source**

Human cells

**Uniprot ID**
[P07141](#)
**Entrez**
[12977](#)
**Expiration Date**

6 months from date of receipt.