
Product Datasheet

Mouse THPO protein (orb707292)

Description

Recombinant Mouse Thrombopoietin Protein

and the target g...

Reactivity

Mouse

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

THPO

Preservatives

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

 Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

 HHHHHHSPVA PACDPRLLNK LLRDSHLLHS
 RLSQCPDVPD LSIPVLLPAV DFSLGEWKTQ
 TEQSKAQDIL GAVSLLLEGV MAARGQLEPS
 CLSLLGQLS GQVRLLLGAL QGLLGTQLPL
 QGRRTAHKDP NALFSLQQL LRGKVRFLLL
 VEGPTLCVRR TLPTTAVPSS TSQLLTLNKF
 PNRTSGLLET NFSVTARTAG PGLLSRLQGF
 RVKITPGQLN QTSRSPVQIS GYLNRTGHPV
 NGTHGLFAGT SLQTLASDI SPGAFNKGSL
 AFNLQGG LPP SPSLAPDGHT PFPSPALPT
 THGSPQLHP LFPDPSTTMP NSTAPHPVTM
 YPHPRNLSQE T

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells