

Product Datasheet

Rat KITLG protein (orb707325)

Description

Recombinant Rat SCF Protein is produced by

encoding Gln...

Reactivity

Rat

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications
Target

KITLG

Preservatives

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MQEICRNPVT DNVKDITKLV ANLPNDYMIT
LNYVAGMDVL PSHCWLRDMV THLSVSLTTL
LDKFSNISEG LSNYSIIDKL GKIVDDLVA
MEENAPKNVK ESLKKPETRN FTPEEFSIF
NRSIDAFKDF MVASDTSDCV LSSTLGPEKD
SRVSVTKPFM LPPVAHHHHH H

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

E. coli

Uniprot ID
[P21581](#)
Entrez
[60427](#)
Expiration Date

6 months from date of receipt.