


Product Datasheet

Alpha-Ketoglutarate Dehydrogenase Microplate Assay Kit (orb707329)

| | |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description | Detection and Quantification of Alpha-Ketoglutarate Dehydrogenase Activity. |
| Reactivity | All |
| Range | 4 µmol/L - 400 µmol/L |
| Tested Applications | FA |
| Storage | Shipped and store at 4°C for 6 months. |
| Note | For research use only |
| Application notes | Detection Method: Colorimetric, Sample Types: Tissue extracts, Cell lysate, Cell culture media, Other biological fluids, Platform: Microplate reader 340 nm |
| Sample Types | Tissue extracts, Cell lysate, Cell culture media, Other biological fluids |
| CAS Number | 9031-02-1 |
| Expiration Date | Please enquire. |



Representative Standard Curve