


Product Datasheet

Resistant Starch Microplate Assay Kit (orb707331)

Description	Detection and Quantification of Resistant Starch Content.
Reactivity	All
Range	0.04 mg/ml - 4 mg/ml
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Tissue extracts, Food, Other agricultural products, Platform: Microplate reader 505 nm
Sample Types	Tissue extracts, Food, Other agricultural products
Expiration Date	Please enquire.



Representative
Standard
Curve