


## Product Datasheet

# Glycogen Branching Enzyme Microplate Assay Kit (orb707332)

<b>Description</b>	Detection and Quantification of Glycogen Branching Enzyme Activity.
<b>Reactivity</b>	All
<b>Range</b>	0.05 mmol/L - 5 mmol/L
<b>Tested Applications</b>	FA
<b>Storage</b>	Shipped and store at 4°C for 6 months.
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Method: Colorimetric, Sample Types: Tissue extracts, Cell lysate, Cell culture media, Other biological fluids, Platform: Microplate reader 660 nm
<b>Sample Types</b>	Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
<b>CAS Number</b>	9001-97-2
<b>Expiration Date</b>	Please enquire.



Representative Standard Curve