


Product Datasheet

Cysteine Microplate Assay Kit (orb707333)

| | |
|----------------------------|---|
| Description | Detection and Quantification of Cysteine Content. |
| Reactivity | All |
| Range | 20 mmol/L - 2000 mmol/L |
| Tested Applications | FA |
| Storage | Shipped and store at 4°C for 6 months. |
| Note | For research use only |
| Application notes | Detection Method: Colorimetric, Sample Types: Tissue extracts, Cell lysate, Cell culture media, Other biological fluids, Platform: Microplate reader 600 nm |
| Sample Types | Tissue extracts, Cell lysate, Cell culture media, Other biological fluids |
| CAS Number | 52-90-4 |
| Expiration Date | Please enquire. |



Representative Standard Curve