


Product Datasheet

NAD-Malate Dehydrogenase Microplate Assay Kit (orb707338)

Description	Detection and Quantification of NAD-Malate Dehydrogenase Activity.
Reactivity	All
Range	4 µmol/L - 400 µmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids, Platform: Microplate reader 340 nm
Sample Types	Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	9001-64-3
Expiration Date	Please enquire.



Representative Standard Curve