


Product Datasheet

Glutathione Peroxidase Microplate Assay Kit (orb707340)

Description	Detection and Quantification of Glutathione Peroxidase Activity.
Reactivity	All
Range	0.01 mmol/L - 5 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Serum, Plasma, Other biological fluids, Tissue extracts, Cell lysate, Cell culture media, Platform: Microplate reader 423 nm
Sample Types	Serum, Plasma, Other biological fluids, Tissue extracts, Cell lysate, Cell culture media
CAS Number	9013-66-5
Expiration Date	Please enquire.



Representative Standard Curve